
**THUASNE**

# Genu Ligaflex®

<b>fr</b>	Genouillère ligamentaire articulée à montants latéraux.....	8
<b>en</b>	Hinged ligament knee brace with lateral uprights .....	10
<b>de</b>	Knieführungsorthese mit Gelenk und seitlichen Schienen .....	12
<b>nl</b>	Scharnierende ligamentaire kniebrace met laterale baleinen.....	15
<b>it</b>	Ginocchiera per legamenti articolata con montanti laterali .....	18
<b>es</b>	Rodillera ligamentaria articulada con montantes laterales.....	20
<b>pt</b>	Joelheira ligamentar articulada com armações laterais.....	23
<b>da</b>	Artikuleret knæbind til ledbånd med stivere i siden.....	25
<b>fi</b>	Nivelletty polven nivelsidetuki, sivupystytuet .....	28
<b>sv</b>	Ledat knästöd för ligamentskada med laterala skenor.....	30
<b>el</b>	Επιγονατίδα συνδέσμων αρθρωτή με πλευρικές ενισχύσεις.....	32
<b>cs</b>	Kolenní vazová ortéza s kloubem a postranními dlahami.....	35
<b>pl</b>	Przegubowa orteza stawu kolanowego ze słupkami bocznymi .....	37
<b>lv</b>	Ceļa locītavas saišu ortoze ar eņģu sānu stiprinājumiem.....	40
<b>lt</b>	Lankstus raiščių antkelis su šoninėmis juostomis .....	42
<b>et</b>	Liigendiga külgtoestusega põlvvetugi .....	45
<b>sl</b>	Pregibna opornica za koleno s stranskimi oporami.....	47
<b>sk</b>	Ohybná ortéza na väzy kolenného kĺbu s bočnými dlahami.....	49
<b>hu</b>	Ízületes térdrögzítő oldalsínekkal .....	52
<b>bg</b>	Съчленена наколенка за връзките на коляното със странични опори.....	54
<b>ro</b>	Orteză de genunchi ligamentară cu articulație și atele laterale.....	57
<b>ru</b>	Лигаментарный ортез для коленного сустава с шарнирами и боковыми креплениями .....	59
<b>hr</b>	Zglobna ortoza za ligamente koljena s bočnim držačima.....	62
<b>zh</b>	帶有橫向支架的關節韌帶護膝.....	64
<b>ar</b>	داعم مفصلي للركبة وأربطتها مزود بدعامتين جانبيتين.....	69

1	32 - 34 cm	37,5 - 43,5 cm
2	35 - 37 cm	40,5 - 46,5 cm
3	38 - 41 cm	44 - 50 cm
3+	38 - 41 cm	48 - 54 cm
4	42 - 45 cm	48 - 54 cm
4+	42 - 45 cm	52 - 58 cm
5	46 - 49 cm	52 - 58 cm
5+	46 - 49 cm	56 - 62 cm
6	50 - 53 cm	56 - 62 cm
6+	50 - 53 cm	60 - 66 cm

1	32 - 34 cm
2	35 - 37 cm
3	38 - 41 cm
4	42 - 45 cm
5	46 - 49 cm
6	50 - 53 cm

<b>fr</b>	Version ouverte
<b>en</b>	Open version
<b>de</b>	Offene Version
<b>nl</b>	Open versie
<b>it</b>	Versione aperta
<b>es</b>	Versión abierta
<b>pt</b>	Versão aberta
<b>da</b>	Åben udgave
<b>fi</b>	Avoin versio
<b>sv</b>	Öppen version
<b>el</b>	Ανοικτή εκδοχή
<b>cs</b>	Otevřená verze
<b>pl</b>	Wersja otwarta
<b>lv</b>	Atvērtā versija
<b>lt</b>	Atviras variantas
<b>et</b>	Avatud tootemudel
<b>sl</b>	Odperto različica
<b>sk</b>	Otvorená verzia
<b>hu</b>	Nyitott változat
<b>bg</b>	Отворена версия
<b>ro</b>	Versiune deschisă
<b>ru</b>	Открытый вариант
<b>hr</b>	Otvorena verzija
<b>zh</b>	开放版本
<b>ar</b>	الاصيغة المفتوحة

<b>fr</b>	Version fermée
<b>en</b>	Closed version
<b>de</b>	Geschlossene Version
<b>nl</b>	Gesloten versie
<b>it</b>	Versione chiusa
<b>es</b>	Versión cerrada
<b>pt</b>	Versão fechada
<b>da</b>	Lukket udgave
<b>fi</b>	Suljettu versio
<b>sv</b>	Stängd version
<b>el</b>	Κλειστή εκδοχή
<b>cs</b>	Uzavřená verze
<b>pl</b>	Wersja zamknięta
<b>lv</b>	Slēgtā versija
<b>lt</b>	Uždaras variantas
<b>et</b>	Kinnine tootemudel
<b>sl</b>	Zaprta različica
<b>sk</b>	Zatvorená verzia
<b>hu</b>	Zárt változat
<b>bg</b>	Затворена версия
<b>ro</b>	Versiune închisă
<b>ru</b>	Закрытый вариант
<b>hr</b>	Zatvorena verzija
<b>zh</b>	封闭版本
<b>ar</b>	الاصيغة المغلقة

<b>fr</b>	Circonférence de la cuisse (14 cm au dessus du genou)
<b>en</b>	Thigh circumference (14 cm above the knee)
<b>de</b>	Oberschenkelumfang (14 cm über dem Knie)
<b>nl</b>	Omvang van de dij (14 cm boven de knie)
<b>it</b>	Circonferenza della coscia (14 cm al di sopra del ginocchio)
<b>es</b>	Circonféncia del muslo (14 cm por encima de la rodilla)
<b>pt</b>	Circonféncia da coxa (14 cm acima do joelho)
<b>da</b>	Lårets omkreds (14 cm over knæet)
<b>fi</b>	Reiden ympärys (14 cm polven yläpuolelta)
<b>sv</b>	Omkrets runt lår (14 cm ovanför knäet)
<b>el</b>	Περίμετρος του μηρού (14 εκ. πάνω από το γόνατο)
<b>cs</b>	Obvod stehna (14 cm nad kolenem)
<b>pl</b>	Obwód uda (14 cm powyżej kolana)
<b>lv</b>	Augšstilba apkārtmērs (14 cm virs ceļgala)
<b>lt</b>	Šlaunies apimtis (14 cm virš kelio)
<b>et</b>	Reie ümbermõõt (14 cm põlvest kõrgemal)
<b>sl</b>	Obseg stegna (14 cm nad kolenom)
<b>sk</b>	Obvod stehna (14 cm nad kolenom)
<b>hu</b>	Comb körfogata (14 cm a térd felett)
<b>bg</b>	Обиколка на бедрото (14 cm над коляното)
<b>ro</b>	Circumferința coapsei (la 14 cm deasupra genunchiului)
<b>ru</b>	Диаметр окружности бедра (на 14 см выше колена)
<b>hr</b>	Opseg bedra (14 cm iznad koljena)
<b>zh</b>	大腿周长 (膝盖上方14厘米处)
<b>ar</b>	محيط الفخذ (14 سم أعلى الركبة)

<b>fr</b>	Circonférence du genou
<b>en</b>	Knee circumference
<b>de</b>	Kniefumfang
<b>nl</b>	Omvang van de knie
<b>it</b>	Circonferenza del ginocchio
<b>es</b>	Circonféncia de la rodilla
<b>pt</b>	Circonféncia do joelho
<b>da</b>	Omkreds af knæ
<b>fi</b>	Polven ympärys
<b>sv</b>	Omkrets runt knä
<b>el</b>	Περίμετρος του γόνατος
<b>cs</b>	Obvod kolena
<b>pl</b>	Obwód kolana
<b>lv</b>	Ceļgala apkārtmērs
<b>lt</b>	Kelio apimtis
<b>et</b>	Põlve ümbermõõt
<b>sl</b>	Obseg kolena
<b>sk</b>	Obvod kolena
<b>hu</b>	Térd körfogata
<b>bg</b>	Обиколка на коляното
<b>ro</b>	Circumferința genunchiului
<b>ru</b>	Диаметр окружности колена
<b>hr</b>	Opseg koljena
<b>zh</b>	膝关节周长
<b>ar</b>	محيط الركبة

<b>fr</b>	Stabilisation	Effet proprioceptif
<b>en</b>	Stabilisation	Proprioceptive effect
<b>de</b>	Stabilisierung	Propriozeptive Wirkung
<b>nl</b>	Stabilisatie	Proprioceptieve werking
<b>it</b>	Stabilizzazione	Effetto proprioceptivo
<b>es</b>	Estabilización	Efecto propioceptivo
<b>pt</b>	Estabilização	Efeito proprioceptivo
<b>da</b>	Stabilisering	Proprioceptisk effekt
<b>fi</b>	Stabiloi	Liikeaistiin kohdistuva vaikutus
<b>sv</b>	Stabilisering	Proprioceptiv effekt
<b>el</b>	Σταθεροποίηση	Αποτέλεσμα ιδιοδεκτικότητας
<b>cs</b>	Stabilizace	Proprioceptivní účinek
<b>pl</b>	Stabilizacja	Propriocepcja
<b>lv</b>	Stabilizēšana	Proprioceptīvs efekts
<b>lt</b>	Stabilizavimas	Proprioceptinis poveikis
<b>et</b>	Stabiliseerimine	Proprioceptiivne toime
<b>sl</b>	Stabiliziranje	Proprioceptivni učinek
<b>sk</b>	Stabilizácia	Proprioceptívny účinok
<b>hu</b>	Stabilizálás	Proprioceptív hatás
<b>bg</b>	Стабилизация	Проприоцептивен ефект
<b>ro</b>	Stabilizare	Efect proprioceptiv
<b>ru</b>	Стабилизация	Проприоцептивный эффект
<b>hr</b>	Stabilizacija	Učinak proprioceptivnog vježbanja
<b>zh</b>		本体感
<b>ar</b>	الاستقرار	مفعول استقبال الحس العميق

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## en HINGED LIGAMENT KNEE BRACE WITH LATERAL UPRIGHTS

### Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Bilateral product.

The product is available in 4 forms: 2 versions (open and closed), each available in 2 lengths (short 30 cm and long 40 cm) to meet the needs of patients (example of open long version in figure A and closed short version in figure B).

Available in 10 sizes for the open version, 6 sizes for the closed version and a version made to measure patient.

The device is composed of:

- an anatomically-shaped knitting integrating a light and soft knit at the back of the knee,
- 2 rigid hinged uprights providing knee stability ①
- the TMS hinge mimicking the natural movements of the knee ②
- 2 anterior half-straps ③, ④ and 2 posterior half-straps ⑤, ⑥ (including a suspension strap ⑦),
- clip buckle system ⑧ for the open version,
- a patellar insert to support patella in the closed version ④.

### Composition

Textile components: polyamide - polyester - polyurethane - elastane - acrylic.

Rigid components: aluminium - polyethylene - polyoxymethylene (open version) - brass (open version) - silicone.

### Properties/Mode of action

Stabilization of knee joint ligaments provided by the rigid hinged side uprights.

An anatomically-shaped knitting for even compression (closed version).

The brace is maintained on the leg thanks to:

- the Suspension Strap ⑧,
  - the silicone-coated threads at the top of the knee brace.
- Staggered straps to prevent tourniquet effect.

### Indications

Conservative treatment of knee ligament injuries.

Joint instability/laxity.

Knee pain and/or swelling management.

Sport activities resumption.

### Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin.

History of venous or lymphatic disorders.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

### Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart; not applicable to the medical devices made to measure patient, the measures are individually taken.

It is recommended that a healthcare professional supervises the first application. Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.

Do not wear the product in a medical imaging device.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

### Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

### Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated.

The straps are numbered using a perforation system: strap ⑧ = 1 perforation, strap ⑨ = 2 perforations... (see diagram for positions).

To avoid damage to the knit fabric, please make sure that the self-gripping pads are always fastened to the designated self-gripping surfaces when putting it on and taking it off.

### For the open version:

- To position the knee brace properly, place the fine knit ① at the back of the knee, with your leg straight; hold the uprights on either side at the hinge and perform a

flexion/extension movement.

- Cover the leg with the knee brace by fastening the touch fasteners at the top and bottom.
- First tighten and fasten the suspension strap ④ located above the calf.
- Clip the buckle of strap ③, and then tighten the strap while holding the buckle lightly.
- Clip and/or tighten straps ③ and ④ successively.

#### **For the closed version:**

- Slip the knee brace on the leg with your foot extended, making sure the patella is properly positioned at the centre of the patellar insert.
- To make it easier to put on, use the pull tabs at the top of the brace and turn it by approx. 90° beforehand.
- First tighten and fasten the suspension strap ④ located above the calf.
- Tighten strap ③, the lower strap of the knee brace.
- Tighten strap ③, the back strap over the thigh.
- Tighten strap ③, the front strap over the thigh.

Tightening the front/back straps adjusts the position of the uprights on the leg to optimize comfort.

If you feel the need to tighten the straps during an activity, make sure you tighten them following the steps indicated in the fitting instructions.

To make the knee brace more comfortable and improve its hold, the uprights can be shaped.

If shaped, consult a healthcare professional.

Apply a slight pressure in the desired direction.

Repeat the process on the other upright if needed.

#### **Care/Maintenance**

Close the self-fastening tabs before washing.

Remove lateral uprights before washing.

To do this:

- Remove the straps, making a note of their original positions.
  - Pull on the top of the product to disengage the upper end of the lateral uprights from of their casings.
  - Remove the uprights via the top opening.
- Replace after washing:
- Insert the hinge via the top opening, with the domed part outwards and flexed towards the back.
  - Insert the upper and lower ends of the uprights into their casings.
  - Check that the loops protrude properly through the openings.
  - Replace the straps, according to their position shown on Fig. A for long version and Fig. B for short version.

Machine washable at 30 °C (delicate programme). If possible use a washing net. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

#### **Storage**

Store at room temperature, preferably in the original packaging.

#### **Disposal**

Dispose of in accordance with local regulations.

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Οι μάντες είναι αφαιρημένοι με ένα σύστημα σφράγισης: μάντας @ = 1 σφραγίδα, μάντας © = 2 σφραγίδες... (βλ. σχέδιο για τις θέσεις).