THUASNE

VENO Elegance

- de Gebrauchsanweisung - Medizinische Kompressionsarmstrümpfe
en Instructions for use - Medical compression sleeves
fr Consignes d'utilisation - Manchons de compression médicale
it Istruzioni per l'uso - Bracciali compressivi medicali
es Instrucciones de utilización - Mangas de compresión terapéuticas para brazos
pt Instruções de uso - Meias de compressão medicinais para braços
nl Gebruiksaanwijzing - Medische armkousen
da Brugsanvisning - Medicinske kompressionsarmstrømper
sv Bruksanvisning - Medicinska kompressionsstrumpor arm
fi Käyttöohje - Lääkinnälliset kompressiohihat
no Bruksanvisning - Medisinske kompresjonsermer
el Οδηγίες χρήσης - Κάλτσες βραχίονα ιατρικής συμπίεσης
cs Návod k použití - Zdravotnické kompresivní pažní návleky
pl Sposób użytkowania - Medyczne rękawy kompresyjne
et Kasutusjuhend - Meditsiinilised kompressioonvarrukad
lv Lietošanas instrukcija - Medicīniskās kompresijas piedurknes
lt Naudojimo instrukcija - Medicininės kompresinės rankovės
sk Návod na použitie - Zdravotnícke kompresívne návleky na paže
sl Navodila za uporabo - Medicinski kompresijski rokav
hu Használati utasítás - Gyógyászati kompressziós karharisnya
ro Instrucțiuni de utilizare - Manșete compresive medicinale pentru brațe
bg Инструкция за употреба - Медицински компресивен чорап за ръка
hr Upute za uporabu - Medicinski kompresijski rukavi
ru Инструкция по использованию - Медицинский компрессионный рукав
vi Hướng dẫn sử dụng - Găng Áp lực Y tế
ko 독일어 - 의료형 팔 압박 스타킹
ja 使用説明書 - 医療用弱圧弾性ストッキング
zh 使用說明 - 医用压力袖套
ar الجوارب الضاغطة الطبية للزراع - جوارب ضغط طبية

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Purpose

Compression sleeves are the basic therapy in the treatment of venous disorders and oedemas. In addition to their medical effect, compression sleeves support your well-being as a whole.

Compression therapy comprises measures of exterior, planar pressure applications in case of venous disorders and lymphatic drainage problems. The application of planar pressure is to prevent the development of oedemas and support venous reflux or lymphatic drainage.

Compression therapy generally is a long-term treatment and requires the corresponding products which – if the medical indication is given – the patient him/herself applies upon the doctor's instruction.

Indications

- Soft indentable lymphoedema (stage I)
- Arm oedema following breast cancer surgery and radiotherapy
- Post-surgery and post-traumatic oedema

After burns and for scar treatment, we recommend our special products for scar therapy.

Absolute contraindications

Peripheral arterial occlusive disease (PAOD) stage III/IV; decompensated heart failure; septic phlebitis; phlegmasia coerulea dolens; acute erysipelas; acute bacterial, viral or allergic inflammations; untreated or advanced dermatoses /contact dermatitis; untreated open wounds; swelling of extremities of unknown origin; untreated malignant tumours.

Relative contraindications

Pronounced exudative dermatoses; malignant lymphoedema; peripheral arterial occlusive disease (PAOD) stage I/II; peripheral neuropathy associated with diabetes mellitus; sensitivity disorders of extremities; compensated heart failure; skin infections; dermatoses (e.g. blistering dermatoses); intolerance of pressure or of product ingredients; primary chronic polyarthritis. In case of open wounds in the treatment area, the wound area must first covered with an adequate wound dressing or a bandage before the compression stockings are put on.

Thuasne Deutschland GmbH accepts no liability for failure to observe these contraindications.

Side effects

No side effects are known if properly used. In case of any negative changes during the application, please consult your doctor or the medical retailer immediately. In case of complaints in connection with the product, such as damage of the knitted fabric or lack of fit, please contact directly your medical retailer.

The manufacturer shall not be liable for damages / injuries which resulted due to improper handling or misuse.

Precautions

Follow the instructions given by your doctor, by an orthopaedic technician or medical expert personnel. If you experience severe pain or discomfort, please see your doctor, an orthopaedic technician or medical expert personnel. Store at room temperature, preferably in the original packaging.

To guarantee its effectiveness and for reasons of hygiene, the product may not be passed on to other patients. Any serious incident that occurs in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient are/is domiciled.

Directions for putting the sleeve on

Put on your compression sleeve in the morning after getting up. Sharp finger nails or jewellery may damage the sleeve. It would be best to put on your compression sleeve with rubber gloves; you can thus roll up the sleeve with the palm of your hand until it fits perfectly without wrinkles. When you use sleeves with an adhesive band, please be sure not to get your hands into the seam of the adhesive band when putting the sleeve on or taking it off.

We recommend our fitting aids if there are difficulties putting on the sleeve.

Care

In order to maintain its function, wash your sleeve every day after its removal. Use a mild detergent for washing by hand or for machine wash (gentle cycle at 40° C) without optical

brightener and fabric softener. In a laundry bag, compression sleeves are specially protected. In order to reduce the drying period, put the compression sleeve into a terry cloth towel, roll it up firmly and press the sleeve out (do not wring it out). Do not leave the moist sleeve in the wet towel. Dry your sleeve in the low cycle of your dryer or simply let it air dry. Never put it in the sun or onto a radiator to dry.

Do not expose the sleeve to any damaging effects, such as: strong sunlight, intense heat, ointments, oils and greases, gasoline and aggressive detergents.

Material composition

Please find the exact information on the sewn-in fabric label.

Note: The sewn-in fabric label is an important component of the compression stockings. If you remove this labelling from a Thuasne product, any claim for warranty, repair or replacement shall be waived.

Storage information and duration of use

Store in a dry place and protect from sunlight. For series products, a maximum shelf life of 36 months applies. This is followed by the duration of use of 6 months for the compression clothing. The hourglass pictogram shows the duration of use printed on the carton label. Customised manufacture concerns products for immediate use.

We guarantee the medical efficacy in series sizes and customised manufacture for a maximum wearing period of six months. This is subject to the correct handling of the sleeves (for example, their care or putting them on or taking them off). If necessary, your doctor should then prescribe a new compression sleeve. In case of any significant change of relevant body sizes (e.g. due to a successful therapy, weight change), new measurements and a new supply may also be justified earlier already.

Disposal

You may dispose of your medical compression sleeves via the residual or household waste. No special criteria exist for their elimination.

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THUASNE (Shanghai) Medical Devices Co., Ltd. · 159 Tianzhou Lu · Shanghai, 200233
Thuasne UK Ltd
THUASNE DEUTSCHLAND GmbH