
**THUASNE**

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**THUASNE**  
120, rue Marius AUFAN  
92300 Levallois-Perret  
France

[www.thuasne.com](http://www.thuasne.com)

**UK Responsible Person (UKRP):**  
Thuasne UK Ltd  
Unit 4 Orchard Business Centre  
North Farm Road  
Tunbridge Wells TN2 3XF,  
United Kingdom

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**en**
**SHOULDER IMMOBILISATION SLING**
**Description/purpose**
This device is a shoulder sling allowing the immobilisation of the shoulder and/or the elbow in internal rotation.
The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.
**Composition**
polyester - polyurethane - polyamide - cotton.
**Properties/Mode of action**
Arm sling with two straps to support the arm and the hand.
Large removable strap to immobilise the arm against the torso.
Bilateral model.
**Indications**
Shoulder and/or elbow injuries or inflammatory conditions.
Conservative treatment for non-surgical proximal humeral fractures.
Shoulder dislocation.
Non surgical fractures of the clavicaula and/or scapula.
Post-operative immobilisation (such as after rotator cuff surgery or dislocation reduction...).
Shoulder subluxation or dislocation prevention (joint instability).
Support for paralytic upper limb.
**Contraindications**
Do not use the product if the diagnosis has not been confirmed.
Do not apply the product in direct contact with broken skin.
Do not use in the event of known allergy to any of the components.
Do not use in certain post-operative conditions (multiple fractures, pin insertion, open fractures, etc.).
Do not use in the event of complex or surgical shoulder fractures, humeral fractures or rib cage fractures.
**Precautions**
Verify the product's integrity before every use.
Do not use the product if it is damaged.
Choose the appropriate size to fit the patient, referring to the size chart.
It is recommended that a healthcare professional supervises the first application.
Strictly comply with your healthcare professional's prescription and recommendations for use.
In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.
For hygiene, security and performance reasons, do not re-use the product for another patient.
Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).
Do not wear the product when driving a vehicle.
Prolonged immobilisation can sometimes cause muscle weakness.
If the support is worn for prolonged periods of time, ensure that blood circulation is not impeded.
It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.
Wear the device above thin clothing.
**Undesirable side-effects**
This device can cause skin reactions (redness, itching, burns, blisters,etc.) or wounds of various degrees of severity.
Any serious incidents related to

the product should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

**Instructions for use/Application**
The help of another person is required to put on this product.
Place the injured forearm in the arm-sling.
Pass the 2 straps over the shoulder.
Cross the 2 straps on the back.
Slide the 2 straps under the injured arm and attach onto the front of the arm-sling.
Attach the large removable strap onto the arm-sling, using the self-fastening closure.
Wrap around the torso and injured arm and attach onto the arm-sling.
**Care**
Product can be washed in accordance with the instructions shown on this leaflet and on the label.
Close the self-fastening tabs before washing.
Machine washable at 30 °C (delicate programme).
The device has been tested over 20 wash cycles.
Do not use detergents, fabric softeners or aggressive products (products containing chlorine).
Do not dry clean.
Do not tumble-dry.
Do not iron.
Squeeze out excess water.
Dry flat.
Dry away from any direct heat source (radiator, sun, etc.).
**Storage**
Store at room temperature, preferably in the original packaging.
**Disposal**
Dispose of in accordance with local regulations.
Keep these instructions for use.

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© Thuasne